


# Bottle Top Dispenser



Bottle Top Dispenser



Bottle Top Dispenser with re-circulation valve



Bottle Top Dispenser with Dual Inlet

Fisherbrand™ offers three lines of bottle top dispensers which caters to every laboratory need from simple dispensing to handling of aggressive and hazardous chemicals and reagents. Fisherbrand dispensers are robust in design with a PTFE piston for smooth movement, easy volume setting, high chemical compatibility and an adjustable delivery nozzle for ease of working in challenging laboratory conditions. These dispensers are fully autoclavable and have easy user re-calibration.

Bottle Top dispenser with re-circulation valve facilitates quick and definite purging (removal of air) without losing any reagent.

Bottle Top Dispenser with Dual Inlet facilitates multiple dispensing modes including rinsing, dilution and bottle re-filling without dismounting the instrument from the reagent bottle.

Whatever your application may be, you can always select a suitable Fisherbrand Bottle Top Dispenser.

Features	Bottle Top Dispenser	Bottle Top Dispenser with re-circulation valve	Bottle Top Dispenser with dual inlet
Dual Inlet Technology			✓
Re-circulation Valve		✓	✓
PTFE Piston	✓	✓	✓
180° Turn Volume Adjustment Knob	✓	✓	✓
Adjustable Delivery Nozzle	✓	✓	✓
360° Rotation of Dispenser on reagent bottle		✓	✓
Individual calibration report with each unit in accordance with ISO 8655 standards.	✓	✓	✓
Five Bottle Adapters (28,32,36,40,45mm)	✓	✓	✓

## Specifications & Ordering Information

Cat. No.	Vol. Range	Increment	Accuracy		CV	
			±%	± ml	±%	± ml
<b>Bottle Top Dispenser</b>						
13-615-10	0.25-2.5 ml	0.05 ml	0.6	0.015	0.2	0.005
13-615-11	0.5-5 ml	0.1 ml	0.6	0.030	0.2	0.010
13-615-12	1-10 ml	0.2 ml	0.6	0.060	0.2	0.020
13-615-13	2.5-30 ml	0.5 ml	0.6	0.180	0.2	0.060
13-615-14	5-60 ml	1.0 ml	0.6	0.360	0.2	0.120
13-615-15	10-100 ml	2.0 ml	0.6	0.600	0.2	0.200
<b>Bottle Top Dispenser with Re-circulation Valve</b>						
13-681-524	0.25-2.5 ml	0.05 ml	0.6	0.015	0.2	0.005
13-681-525	0.5-5 ml	0.1 ml	0.6	0.030	0.2	0.010
13-681-526	1-10 ml	0.2 ml	0.6	0.060	0.2	0.020
13-681-527	2.5-30 ml	0.5 ml	0.6	0.180	0.2	0.060
13-681-528	5-60 ml	1.0 ml	0.6	0.360	0.2	0.120
13-681-529	10-100 ml	2.0 ml	0.6	0.600	0.2	0.200
<b>Bottle Top Dispenser with Dual Inlet</b>						
13-681-530	0.25-2.5 ml	0.05 ml	0.6	0.015	0.2	0.005
13-681-531	0.5-5 ml	0.1 ml	0.6	0.030	0.2	0.010
13-681-532	1-10 ml	0.2 ml	0.6	0.060	0.2	0.020
13-681-533	2.5-30 ml	0.5 ml	0.6	0.180	0.2	0.060
13-681-534	5-60 ml	1.0 ml	0.6	0.360	0.2	0.120
13-681-535	10-100 ml	2.0 ml	0.6	0.600	0.2	0.200

### In the United States:

For customer service, call 1-800-766-7000  
 To fax an order, use 1-800-926-1166  
 To order online: fishersci.com

### In Canada:

For customer service, call 1-800-234-7437  
 To fax an order, use 1-800-463-2996  
 To order online: fishersci.ca

